3MSM Health Care Academy

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
Aerobic Bacteria Count: 135

Count all red colonies regardless of size or intensity.

AC
Aerobic Count Plate

Product Code	6400 (box of 100 units) 6406 (box of 1,000 units)
Method: AOAC® Official Method of Analysis™	Milk and Dairy Products 48h ± 3h at 32°C ± 1°C All Other Foods 48h ± 3h at 35°C ± 1°C
NF Validation by AFNOR Certification	For All Foods 72h ± 3h at 30°C ± 1°C All Foods Except Dairy and Raw Shellfish 48h ± 3h at 30°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.6–7.2
Recommended Counting Range	<300 CFU
Counting Area	20cm²
Spreader Type	3M™ Petrifilm™ Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

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
Aerobic Bacteria Count: 88

Blue and red indicator dyes in the plate color the colonies. Count all colonies regardless of their size or color intensity.

Product Code	6478 (box of 50 units) 6479 (case of 500 units)
Method: AOAC® Official Method of Analysis™	Dairy (Not Including Powders) and Seafood 24h ± 2h at 32°C ± 1°C All Other Foods 24h ± 2h at 35°C ± 1°C Dairy Powders Including Whey Powders 48h ± 3h at 32°C ± 1°C
NF Validation by AFNOR Certification	Non-Powdered Dairy Products 28h ± 2h at 30°C ± 1°C Powdered Milk 48h ± 3h at 30°C ± 1°C
Incubation	Stacks of 40 or less (SMEDP: 20 or less)
Optimum pH Range	>5.0
Recommended Counting Range	<300 CFU
Counting Area	30cm²
Spreader Type	3M™ Petrifilm™ Flat Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

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
Enterobacteriaceae Count: 13

Enterobacteriaceae are red colonies with yellow zones and/or red colonies with gas bubbles with or without yellow zones.

EB
Enterobacteriaceae Count Plate

Product Code	6420 (box of 50 units) 6421 (box of 1,000 units)
Method: AOAC® Official Method of Analysis™	Selected Foods 24h ± 2h at 37°C ± 1°C
NF Validation by AFNOR Certification	All Human Food Products, Animal Feed and Industrial Environmental Samples 24h ± 2h at 30°C ± 1°C or 37°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.5–7.5
Recommended Counting Range	<100 CFU
Counting Area	20cm ²
Spreader Type	3M™ Petrifilm™ Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

EB
Enterobacteriaceae Count Plate

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
AOAC® OMA, NF Validation Certified Method 3M 01/02–09/89 B
Coliform Count: 69 (colonies with gas)

NF Validation by AFNOR Certified Method 3M 01/02–09/89 A and 01/02–09/89 C
Coliform Count: 97 (total colonies)

CC**Coliform Count Plate**

Product Code	6410 (box of 50 units) 6416 (box of 1,000 units)
Method: AOAC® Official Method of Analysis™	Milk and Other Dairy Products 24h ± 2h at 32°C ± 1°C Foods 24h ± 2h at 35°C ± 1°C
NF Validation by AFNOR Certification	3M 01/02-09/89 B All Human Food Products (Except Raw Shellfish) 24h ± 2h at 30°C ± 1°C or 37°C ± 1°C
3M 01/02–09/89 A All Human Food Products (Except Shellfish), Pet Food and Environmental Samples 24h ± 2h at 30°C ± 1°C or 37°C ± 1°C	3M 01/02–09/89 C All Human Food Products 24h ± 2h at 44°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.6–7.2
Recommended Counting Range	<150 CFU
Counting Area	20cm ²
Spreader Type	3M™ Petrifilm™ Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

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
Coliform Colony Enumeration
(8–24 hours)

RCC
Rapid Coliform Count Plate

Product Code	6402 (box of 50 units) 6412 (box of 500 units)
Method: AOAC® Official Method of Analysis™	Foods Up to 24h ± 2h at 35°C ± 1°C
NF Validation by AFNOR Certification	Processed Pork and Seafood 14h ± 30min at 30°C ± 1°C or 24h ± 2h at 30°C ± 1°C All Other Human Foods 14h ± 30min at 35°C ± 1°C or 24h ± 2h at 35°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.5–7.5
Recommended Counting Range	<150 CFU
Counting Area	20cm ²
Spreader Type	3M™ Petrifilm™ Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

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
HSCC
High-Sensitivity Coliform Count Plate

Product Code	6405 (box of 50 units) 6415 (box of 500 units)
Method: AOAC® Official Method of Analysis™	Dairy Products 24h ± 2h at 32°C ± 1°C
NF Validation by AFNOR Certification	Human Food Products 24h ± 2h at 30°C ± 1°C or 37°C ± 1°C
Method	Thermotolerant Coliforms 24h ± 2h at 44°C ± 1°C Incubator humidification is required at this elevated temperature.
Incubation	Stacks of 10 or less
Optimum pH Range	6.5–7.5
Recommended Counting Range	<150 CFU
Counting Area	60cm ²
Spreader Type	3M™ Petrifilm™ High Sensitivity Spreader (6481)

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

HSCC
High-Sensitivity Coliform Count Plate

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
E. coli Count: 49 (blue colonies with gas)
Total Coliform Count: 87 (red and blue colonies with gas)

EC
E. coli/Coliform Count Plate

Product Code	6404 (box of 50 units) 6414 (box of 500 units)
Method: AOAC® Official Method of Analysis™	Coliform in Foods 24h ± 2h at 35°C ± 1°C E. coli in Foods 48h ± 4h at 35°C ± 1°C Poultry, Meats and Seafood Coliforms and E. coli 24h ± 2h at 35°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.6–7.2
Recommended Counting Range	<150 CFU
Counting Area	20cm ²
Spreader Type	3M™ Petrifilm™ Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

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

E. coli Count: 56

SEC
Select *E. coli* Count Plate

Product Code	6434 (box of 50 units) 6435 (box of 500 units)
NF Validation by AFNOR Certification	All Human Food Products, Pet Food and Industrial Environmental Samples 24h ± 2h at 42°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.5–7.5
Recommended Counting Range	<150 CFU
Counting Area	20cm ²
Spreader Type	3M™ Petrifilm™ Spreader

Sold in select regions only.

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

SEC
Select *E. coli* Count Plate

3M™ Petrifilm™


- Staph Express

***S. aureus* Count: 116**
Red-violet colonies are *S. aureus*.

STX
Staph Express Count Plate

Staph Express Confirmation Disk

Pink zones form with DNase reaction from *S. aureus*.

© 3M 2022. All Rights Reserved. 3M Confidential.


Product Code	Plates 6490 (box of 50 units) 6491 (box of 500 units) Disks 6492 (box of 20 units) 6493 (box of 100 units)
Method: AOAC® Official Method of Analysis™	In Food Plate: 24h ± 2h at 35°C ± 1°C or 37°C ± 1°C Disk: 1–3h at 35°C ± 1°C or 37°C ± 1°C
NF Validation by AFNOR Certification	All Human Food and Pet Food Plate: 24h ± 2h at 37°C ± 1°C Disk: 3h at 37°C ± 1°C
Incubation	Stacks of 20 or less
Optimum pH Range	6.0–8.0
Recommended Counting Range	<150 CFU
Counting Area	30cm ²
Spreader Type	3M™ Petrifilm™ Flat Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

STX
Staph Express Count Plate

3M™ Petrifilm™

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YM
Yeast and Mold Count Plate


Product Code	6407 (box of 100 units) 6417 (box of 1,000 units)
Method: AOAC® Official Method of Analysis™	Foods 5 days at 20–25°C
Incubation	Stacks of 20 or less
Optimum pH Range	Not Applicable
Recommended Counting Range	<150 CFU
Counting Area	30cm ²
Spreader Type	3M™ Petrifilm™ Yeast and Mold Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

YM
Yeast and Mold Count Plate

3M™ Petrifilm™

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Mold Count: 12 (after 48 hours)

Mold Characteristics: Large colonies, colonies have diffused edges, blue-green to variable upon prolonged incubation, colonies appear flat and colonies have a dark center with diffused edge.

Product Code	6475 (box of 50 units) 6477 (case of 500 units)
Method: AOAC® Official Method of Analysis™	Select Foods 48–60h at 25°C ± 1°C or 28°C ± 1°C
NF Validation by AFNOR Certification	All Human Food Products, Animal Feed and Environmental Products (Primary Production Samples Excepted) 60–72h at 25°C ± 1°C or 28°C ± 1°C
Incubation	Stacks of 40 or less
Optimum pH Range	Not Applicable
Recommended Counting Range	<150 CFU
Counting Area	30cm²
Spreader Type	3M™ Petrifilm™ Flat Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

RYM
Rapid Yeast and Mold Count Plate

RYM
Rapid Yeast and Mold Count Plate

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
Circle 1: Homofermentative colony (non-gas producing).

Circle 2: Heterofermentative colony (gas producing).

LAB
Lactic Acid Bacteria Count Plate

Product Code	6461 (box of 50 units) 6462 (case of 500 units)
Method	Foods and Environmental Samples Incubate 48h ± 3h at appropriate temperature between 28–37°C
Incubation	Stacks of 20 or less
Optimum pH Range	Not Applicable
Recommended Counting Range	Colonies Without Gas <300 CFU Colonies With and Without Gas <150 CFU
Counting Area	30cm ²
Spreader Type	3M™ Petrifilm™ Flat Spreader

Visit the 3M Food Safety website at [3M.com/FoodSafety](https://www.3m.com/FoodSafety) for product information and product instructions which include information regarding validation recognition, scope and end of validity.

LAB
Lactic Acid Bacteria Count Plate

Questions?